

866,778,7826

Global Equipment Supplier to the Pharmaceutical and Supplement Industries BEIJING



High Speed Automatic

Encapsulation



- Pocket 400 & 800
- CFM-1200
- CFM-2000
- CFM-3500
- CFM-7500

Robust Repeatable Reliable

Right Stuff Equipment 1605 E. 69th Ave. Denver, CO 80229



866 778 7826

Global Equipment
Supplier to the
Pharmaceutical and
Supplement
Industries
BEIJING



Pocket 400 & 800



- Reliable and quiet indexed gearbox drive
- Small footprint
- Entry level pricing
- PLC and touch screen control
- Meets cGMP and FDA requirements

Max output (Cap/min)	400	800
No. of Segment Bores	3	7
Suitable Capsule Size	00#~5#	00#~5#
Filling Precision	±3%	±3%
Noise	≤70dB	≤70dB
Main Motor Power	.75 Kw (1HP)	.75 Kw (1HP)
Dimensions (L x W x H)	800 x 970 x 1870 mm (31"x38"x73")	800 x 970 x 1870 mm (31"x38"x73")
Net Weight	850Kg (1900 lbs.)	850Kg (1900 lbs.)

- A pellet attachment is available to fill capsules with pellets, tablets or granules
- Cantilever type HMI touch screen w/ intuitive Siemens controls (AB available)
- Precise dosage with dynamic adjustments
- Meets cGMP and FDA requirements

	CFM-1200	CFM-2000
Max output (Capsules/min)	1200	2000
No. of Segment Bores	9 single line	9 double line
Suitable Capsule Size	00#~5#	00#~5#
Filling Precision	±3%	±3%
Noise	≤70dB	≤70dB
Total Power	5.5Kw (7.3HP)	7.5Kw (10HP)
Dimensions	1200 x 1014 x 2000	1310 x 1280 x 2000
$(\mathbf{L} \mathbf{x} \mathbf{W} \mathbf{x} \mathbf{H})$	(47"x40"x79")	51" x 45" x 79"
Weight	1600 Kg (3500lbs)	1700 Kg (3750 lbs.)

CFM 1200 & CFM 2000



Right Stuff Equipment 1605 E. 69th Ave. Denver, CO 80229



866 778 7826

Global Equipment Supplier to the Pharmaceutical and Supplement Industries BEIJING





- Built with a precise drive train and gear box producing a smooth and reliable operation for a long life.
- Balance of high production output with high filling precision and flexibility
- Robust and repeatable operation
- Standard controls include Seimens PLC and 10"HMI

Max output (Capsules/min)	3500
No. of segment bores	24 double line
Suitable Capsule Size	00#~5#
Filling precision	±3%
Noise	≤70dB
Total Power	7.5Kw (10HP)
Dimensions (L x W x H)	1310 x 1280 x 2000 51" x 50" x 79"
Net Weight	1850 Kg (4078 lbs.)

Right Stuff Equipment 1605 E. 69th Ave. Denver, CO 80229



866 778 7826

Global Equipment Supplier to the Pharmaceutical and Supplement Industries BEIJING





- Excellent cost to performance ratio
- Flexible machine set up
- Ultra-high production capacity reduces labor costs
- Suitable for mass production of a single product
- Control unit and capsule filler are separated allowing easier placement in the production room and reduces interference of vibration from CFM

Max output (Capsules/min)	7500
No. of Segment Bores	52/4
Suitable for Capsule Size	00#~5#
Filling Precision	±3%
Noise	≤70dB
Total Power	10kW (13.4hp)
Net Weight	3000 Kg (6613 lbs.)
Dimensions	1738 x 1226 x 2000
$(\mathbf{L} \mathbf{x} \mathbf{W} \mathbf{x} \mathbf{H})$	68" x 48" x 79"

Right Stuff Equipment 1605 E. 69th Ave. Denver, CO 80229



866 778 7826

Global Equipment Supplier to the Pharmaceutical and Supplement Industries BEIJING





D' '	1160 1710 1900	
Dimension	1160 x 1710 x 1800 mm	
	(45" x 67" x 71")	
Speed (tablets/min)	0-7000	
Net Weight	700 Kg (1543 lbs.)	
Filling quantities per	8-9999	
bottle		
Tablets Specifications	Thickness: 3-10mm	
•	Diameter.:3-15mm	
	Length 3-25mm	
Number of Tracks:	2x8	
Bottle Size:	Diameter 30- 75	
	Height: 40-250mm	
Hopper Cubage:	2 x 20L	
Compressed Air:	0.5 Mpa	
Conveyor Length:	2400mm (95")	
Number of Vibrators	6	

- 16 Counting Channels
- 6-step electromagnetic vibrator
- Infrared Dynamic Scan Technology
- Rejection of incorrect count bottles
- Easy to disassemble and clean
- 304 Stainless Steel Construction

GLP 800 Bottle Unscrambler



LP 700A Bottle Unscrambler



Right Stuff Equipment 1605 E. 69th Ave. Denver, CO 80229



866 778 7826

Global Equipment
Supplier to the
Pharmaceutical and
Supplement
Industries
BEIJING







Providing Equipment Sales, Service and Engineering Support to the Food, Beverage, Supplement and Pharmaceutical Industries Denver CO

(866)778-7826

Right Stuff Equipment 1605 E. 69th Ave. Denver, CO 80229